

Product Specification Sheet

**Tubby Like Protein-1 (TULP1) Antibodies**

|                       |   |                     |
|-----------------------|---|---------------------|
| <b>Cat # TULP11-S</b> | Rabbit Anti-Human TULP1 antisera # 1            | <b>Size:</b> 100 ul |
| <b>Cat # TULP11-A</b> | Rabbit Anti-Human TULP1 IgG # 1 (Affinity Pure) | <b>Size:</b> 100 ug |
| <b>Cat # TULP11-P</b> | Human TULP1 Control peptide # 1                 | <b>Size:</b> 100 ug |

Several common diseases such as type II diabetes, hypertension, cardiovascular diseases, hyperlipidemia, and some cancers are associated with obesity. An abnormal increase in body fat relative to lean tissue mass has been used as an indicator of obesity. High fat diet, certain environmental factors, and genetic linkage are the primary causes of obesity. In order to understand the genetic basis of obesity, several monogenic murine obesity models have been characterized including *obese (Ob)*, *diabetes (db)*, *fat (fat)*, *agouti yellow (A<sup>y</sup>)*, and *tubby (tub)*. More recently, *Tub*, the human homolog of mouse *Tub*, **TULP1** & **TULP2** (for **Tubby Like Proteins**) and Agouti related protein (AGRP) have been cloned. The obesity associated with *Ay* mice may be due to ectopic expression of a secreted protein Agouti. *Agouti* protein (132 aa in human) is normally expressed in skin but its ubiquitous expression causes obesity. Agouti is a paracrine-signaling molecule that affects pigmentation by inhibiting the melanocortin receptor 1 (MCR-1). However, recombinant Agouti protein also antagonizes the MC2R and MC4R. AGRP (132 aa in human, chromosome 16q21) is normally expressed in adrenal and hypothalamus. AGRP levels are increased several folds in *ob/ob* mice. AGRP is a strong antagonist of MC3R and MC4R. Ubiquitous expression of AGRP in transgenic mice causes obesity without altering skin pigmentation.

**Source of Antigen and Antibodies**

|                        |  |
|------------------------|--|
| <b>Antigen</b>         | 18aa peptide of Human TULP1 (1); <b>Designated (TULP11-P or control peptide)</b> epitope location ~ Mid-region |
| <b>Ab Host/type</b>    | Rabbit, polyclonal Unpurified antiserum (cat # TULP11-S) Aff pure IgG (cat # TULP11-A)                         |
| <b>2-ab</b>            | <b>Goat Anti-rabbit IgG-HRP</b> cat # 20320 (AP, biotin, FITC conjugates also available)                       |
| <b>-ve control IgG</b> | <b># 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as –ve control</b>           |

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied 0.05% azide, **Reconstitute** powder in 100 ul PBS

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder** in PBS at 1mg/ml

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at -200C and powder at 4oC or -20oC..

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Shipping:** 4oC for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure using Chemiluminescence technique).

**ELISA** (1:10K-1:100K; using 50-100 ng of control peptide/well).

**Histochemistry & Immunofluorescence:** not tested. We recommend the use of affinity pure antibody at 2-20 ug/ml.

**Specificity & Cross-reactivity**

The human TULP11-P peptide has no significant homology with mouse *tubby*, *Agouti*, *AGRP*, and *TULP2*. Antibody crossreactivity in various species is not known. The control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at the web site).

**General References:** (1) North MA et al (1997) PNAS 94, 3128-3133

*\*This product is for In vitro research use only.*

**Related material available from ADI**

Anti-Agouti, AGRP, *Tubby*, *TUB*, *TULP1*, *TULP2*, Leptin, and Melanocortin receptors 91-5)

Anti-Rabbit IgG-HRP Conjugate and ECL Reagents

**Western Blot Recycling Kit (Strips blots in 5 minutes)** and re-use the same blot with multiple antibodies

TULP11-S-A-P

71219A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)