


Product Specification Sheet

Tai Forest virus Glycoprotein (TAFV GP, Cote d'Ivoire-94) Antibodies and Controls

| | | |
|--|---|---------------------|
| <input type="checkbox"/> Cat # TVGP11-S | Rabbit Anti-Tai Forest virus Glycoprotein (TAFV GP, Cote d'Ivoire-94) antiserum | SIZE: 100 ul |
| <input type="checkbox"/> Cat # TVGP11-A | Rabbit Anti-Tai Forest virus Glycoprotein (TAFV GP, Cote d'Ivoire-94) IgG | SIZE: 100 ul |
| <input type="checkbox"/> Cat # TVGP11-C | Recombinant Tai Forest virus Glycoprotein (TAFV GP, Cote d'Ivoire-94) control for western | SIZE: 100 ul |

Ebola virus (EBOV), formerly Zaire ebolavirus causes severe disease in humans and in nonhuman primates in the form of viral hemorrhagic fever. Zaire ebolavirus is a virological taxon included in the genus Ebolavirus, family Filoviridae, order Mononegavirales. The species has a single virus member, Ebola virus (EBOV). **Ebolavirus species Zaire (ZEBOV)** causes highly lethal hemorrhagic fever, resulting in the death of **90%** of patients within days. Most information on immune responses to ZEBOV comes from in vitro studies and animal models. Ebola Zaire attacks every organ and tissue in the human body except skeletal muscle and bone. Ebola is classified as a **Level 4** pathogens (higher than AIDS) with a 2 to 21 day (7 to 14 days' average) incubation period. There are currently four known strains of Ebola: Zaire, Sudan, Reston and Tai. All cause illness in sub-human primates. Only Ebola Reston does not cause illness in humans. The mortality rate of Ebola victims is between 60% and 90%; with Ebola Sudan at 60% and Ebola Zaire at 90%.

The virions are tubular in general form but variable in overall shape and may appear as the classic shepherd's crook or eyebolt. Ebola virions consist of 7 structural proteins. At the center is the helical **ribonucleocapsid**, which consists of the genomic RNA wrapped around a polymer of **nucleoproteins (NP)**. Associated with the ribonucleoprotein is the RNA-dependent **RNA polymerase (L)**



with the **polymerase cofactor (VP35)** and a **transcription activator (VP30)**. The ribonucleoprotein is embedded in a matrix, formed by the major (VP40) and minor (VP24) matrix proteins. They are surrounded by a **lipid membrane** derived from the host cell membrane. The membrane anchors a glycoprotein (GP_{1,2}) that projects 7 to 10 nm spikes away from its surface. While nearly identical to Marburg **virions** in structure, ebola virions are antigenically distinct. The most common diagnostic methods are RT-PCR in conjunction with antigen-capture ELISA which can be performed in field or mobile hospitals and laboratories. There are currently no FDA-approved vaccines for the prevention of EVD. The most promising ones are DNA vaccines or are based on adenoviruses, vesicular **stomatitis Indiana virus (VSIV)** or **filovirus-like particles (VLPs)** as all of these candidates could protect nonhuman primates from Ebola virus-induced disease. DNA vaccines, adenovirus-based vaccines, and VSIV-based vaccines have entered clinical trials.

Source of Antigen and Antibodies

| | |
|------------------------|--|
| Antigen | Recombinant protein TAFV GP (33-650aa; <i>Tai Forest ebolavirus (strain Cote d'Ivoire-94) (TAFV) (Cote d'Ivoire Ebola virus)</i>) |
| Ab Host/type | Rabbit, Polyclonal IgG antiserum (Cat#TVGP11-S) or purified (Cat# TVGP11-A). Supplied in 0.05% Sodium Azide as preservative. |
| 2-Ab | Goat Anti-rabbit IgG-HRP cat # 20320 (AP, biotin, FITC conjugates) |
| -ve control IgG | Cat# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control |

Cat#TVGP11-C, positive control

Tai Forest Ebola Virus glycoprotein was expressed in E. coli cells as his-tag fusion protein (33-650 aa, >95%, ~70 kDa). Purified glycoprotein for Western blot +ve control (Cat# TVGP11-C) is supplied in SDS-PAGE sample buffer (reduced). Load 10 ul/lane of # TVGP11-C for good visibility with antibody Cat # TVGP11-S. Store at -20oC in suitable size aliquots. SDS may crystallize in cold conditions. It should be re dissolved by warming before taking it from the stock. It should be heated once prior to loading on gels. If the product has been stored for several weeks, then it may be preferable to add 5 ul of fresh 2x sample buffer per 10 ul of the

TVGP11-C solution prior to heating and loading on gels. This preparation is not biologically active. It is not suitable for ELISA or other applications where native protein is required. Do not freeze, thaw, or heat repeatedly.

Form & Storage of Antibodies/Peptide Control

Antibody
 100 ul solution lyophilized powder
 Buffer: PBS+0.05% azide **Reconstitute powder** 100 ul water

Storage
Short-term: unopened, undiluted vials for less than a week at 4oC.
Long-term: at -20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

Stability: 6-12 months at -20oC or below.
Shipping: 4oC for solutions and room temp for powder.

Recommended Usage

Western Blotting: Assay-dependent dilution (range 1:200-:1000). The antibody detected recombinant Tai Forest Ebola virus GP when tested at 100-200 ng/ml. ~70 kda recombinant. Use ADI #TVGP11-C as western control.

ELISA: Assay dependent. (1:1000-10,00; using 100-200 ng antigen/well). Use ADI #TVGP15-R-10.

Histochemistry & Immunofluorescence: not tested.

Specificity & Crossreactivity: This antibody reacts with the Tai Forest GP strain of the Ebola virus. Ebola GP are 54-55% related (Zaire Ebola GP protein is significantly conserved in various serotypes: Bundibugyo (65%), Tai Forest D'Ivoire (64%), Sudan (54%), Reston (57%). Antibodies and recombinant GP proteins from various Ebola and Marburg viruses are available for control studies.

References: Thomas W (2010) Archives of Virology 155 (12): 2083-2103. Taylor D (2010) BMC Evolutionary Biology 10: 193. Feldmann H (2005) A. Virus Taxonomy—Eighth Report of the International Committee on Taxonomy of Viruses. 645-653.

*This product is for in vitro research use only.

Related material available from ADI

| Catalog# | Prod Description |
|--------------|--|
| TVGP11-A | Anti-Ebola virus glycoprotein (Recombinant) IgG, purified |
| TVGP11-C | Rec. (sf9) Zaire-Ebola virus glycoprotein protein control WB |
| EVGP15-A | Anti-Zaire Ebola virus glycoprotein (GP, 1-676aa/DNA vaccine) IgG, |
| EVGP16-A | Anti-Zaire Ebola virus glycoprotein (GP 1-652aa/DNA vaccine) IgG |
| EVGP16-R-10 | Rec. (sf9) Sudan-Ebola virus glycoprotein (minus transmembrane domain, his-tag, 68 kDa), purified |
| EVGP17-R-10 | Recombinant (sf9) Zaire-Ebola virus glycoprotein (minus transmembrane domain, his-tag, 68 kDa), purified |
| EVP401-A | Anti-Zaire-Ebola virus VP40 peptide (EBOV VP40) IgG, |
| EVP401-C | Rec. Zaire-Ebola virus VP40 protein control for Western |
| EVP405-R-10 | Rec. (E.coli) Zaire-Ebola virus VP40 |
| AE-320520-1 | Human Anti-Ebola virus Nucleoprotein (NP) IgG ELISA Kit, |
| AE-320530-1 | Human Anti-Ebola virus Nucleoprotein (NP) IgM ELISA Kit |
| AE-320620-1 | Human Anti-Zaire-Ebola virus glycoprotein (GP) IgG ELISA |
| AE-320720-1 | Human Anti-Zaire-Ebola virus VP40 IgG ELISA Kit, |
| AE-320730-1 | Human Anti-Zaire-Ebola virus VP40 IgM ELISA Kit, |
| AE-320800-48 | Human Zaire-Ebola Virus antigen ELISA Kit, 48 tests, |
| AE-320800-96 | Human Zaire-Ebola Virus antigen ELISA Kit, 96 tests, |

TVGP11-A-Recombinant-Ebola-Tai-Forest-GP 151210A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com