

Product Specification Sheet

Mouse Uncoupling Protein 2 (UCP2)

<input type="checkbox"/> Cat. UCP25-R-10	Recombinant purified Mouse UCP2 protein	SIZE: 10 ug
<input type="checkbox"/> Cat. UCP25-R-50	Recombinant purified Mouse UCP2 protein	SIZE: 50 ug

The regulation of body weight depends upon the calorie intake and expenditure. White and brown adipose tissues (BAT and WAT, respectively) play a central role in body weight and energy expenditure. WAT is the major site for energy storage via triglyceride synthesis and mobilization via lipolysis. **Uncoupling proteins (UCP1-5)** are a family of mitochondria transport proteins that play a critical role in thermoregulatory heat production and maintenance of basal metabolic rate. BAT is able to dissipate energy as heat via uncoupled mitochondrial respiration by a mitochondrial anion carrier, uncoupling protein 1 (UCP1). UCP1 is predicted to contain 6 trans-membrane (TM) domains, a putative purine nucleotide-binding domain (PNBD) and three-mitochondrial energy transfer protein domains (ETPDs). Both amino and C-termini are predicted to be cytoplasmic. Mouse/rat UCP1 is A 307 AA mitochondrial uncoupling protein (1). It is only 59% homologous with UCP2 (2). UCP1 is primarily expressed in brown adipose tissues.

Mouse/rat **UCP2** is A 309 AA (human 309 aa chromosome 7; ~95% homology) mitochondrial uncoupling protein (1). It is only 59% homologous with UCP1 found in brown adipose tissues. UCP2 has wide tissue distribution in mouse tissues (brain, kidney, liver, brown adipose tissue, heart, and muscle). UCP2 may play a critical role in energy balance, body weight, and thermoregulation (1, 2).

Source of Antigen and Antibodies

Mouse UCP2 protein (full length; gene accession # NM_011671) was expressed as fusion protein (His tag- UCP2) in E.coli and purified (>95% with major band at ~32kDa). Purified Mouse UCP2 protein is supplied in 50mM Tris-HCl, pH 8.0, 250mM NaCl, 1.5% sodium lauroyl Sarcosine, 1mM EDTA, 10mM beta-mercapto-ethanol at a conc. of 1 mg/ml (lot specific conc. is provided on the vial). Store at -20°C or below. Avoid frequent freeze and thaw.

Recommended Usage

For Western load, 50-200 ng protein per lane and probe with appropriate antibodies.
For ELISA, coat at 0.2-2 ug/ml.

General References: (1) Kozak LP (1988) JBC 263, 12274; Bouillaud F (1986) JBC 261, 1487; Ridley RG (1986) Nucl. Acid Res. 14, 4025; Miroux B (1993); EMBO J 12, 3739; Cassard AM et al (1990) JCB 43, 255; (2) Fleury C (1997) Nature Genetics 15, 269;

(2) Citations of ADI's Antibodies (see web site for updated list)
Mulligan, 2007, J Physiology, 580, 677-684, Apr 2007; Jonathan B, 2007, PNAS; 104: 5680 – 5685; Rodriguez E, 2001, Eur. J. Physiology 442, 396-403; Ryu J-W, 2003, BBRC 303, 726-731; Nibbelink M, 2001, J. Biol. Chem. 276, 47295; Klein J, 2004, J. Endocrinol., 183: 299 – 307; Lee S-H, 2004, J. Lipid Res., 45, 1674-1682; Justo R, 2005, Life Sciences, 76, 1147-1158;
See more refs at the web site.

Related Items

Catalog#	ProdDescription
UCP11-A	Anti-Mouse Uncoupling Protein 1 (UCP1) IgG # 1, aff pure
UCP11-P	Mouse Uncoupling Protein 1 (UCP1) Control/blocking
UCP12-A	Anti-Mouse Uncoupling Protein 1 (UCP1) IgG #2, aff pure
UCP12-P	Mouse Uncoupling Protein 1 (UCP1) Control/blocking
UCP15-R-10	Recombinant purified Mouse Uncoupling Protein 1 (UCP1) (full length, His-tag)
UCP16-R-10	Recombinant purified Human Uncoupling Protein 1 (UCP1) (full length, His-tag)

UCP21-A	Anti-Mouse Uncoupling Protein 2 (UCP2) IgG # 1, aff pure
UCP21-P	Mouse Uncoupling Protein 2 (UCP2) Control/blocking peptide #1
UCP22-A	Anti-Mouse Uncoupling Protein 2 (UCP2) IgG # 2, aff pure
UCP22-P	Mouse Uncoupling Protein 2 (UCP2) Control/blocking
UCP23-A	Anti-Human Uncoupling Protein 2 (UCP2) IgG # 3, aff pure (UCP2) (full length, His-tag) WB +ve control
UCP23-P	Human Uncoupling Protein 2 (UCP2) Control/blocking
UCP25-R-10	Recombinant purified Mouse Uncoupling Protein 2 (UCP2) (full length, His-tag)
UCP26-R-10	Recombinant purified Human Uncoupling Protein 2 (UCP2) (full length, His-tag)

UCP31-A	Anti-Human Uncoupling Protein 3 (UCP3) IgG # 1, aff
UCP31-P	Human Uncoupling Protein 3 (UCP3) Control/blocking
UCP31-S	Anti-Human Uncoupling Protein 3 (UCP3) antiserum # 1
UCP32-A	Anti-Human Uncoupling Protein 3 (UCP3) IgG # 2, aff pure
UCP32-S	Anti-Human Uncoupling Protein 3 (UCP3) antiserum # 2
UCP33-A	Anti-Rat Uncoupling Protein 3 (UCP3) IgG # 3, aff pure
UCP33-P	Rat Uncoupling Protein 3 (UCP3) Control/blocking peptide
UCP34-A	Anti-Human Uncoupling Protein 3 (UCP3) IgG # 4, aff
UCP34-P	Human Uncoupling Protein 3 (UCP3) Control/blocking
UCP35-R-10	Recombinant purified Mouse Uncoupling Protein 3 (UCP3) (full length, His-tag)
UCP36-R-10	Recombinant purified Human Uncoupling Protein 3 (UCP3) (full length, His-tag)

UCP41-A	Anti-Human Uncoupling Protein 4 (UCP4) IgG # 1, aff pure
UCP41-P	Human Uncoupling Protein 4 (UCP4) Control/blocking
UCP42-A	Anti-Rat Uncoupling Protein 4 (UCP4) IgG # 2, aff pure
UCP42-P	Rat Uncoupling Protein 4 (UCP4) Control/blocking peptide
UCP42-S	Anti-Rat Uncoupling Protein 4 (UCP4) antiserum # 2

UCP51-A	Anti-Rat Uncoupling Protein 5 (UCP5) IgG #1, aff pure
UCP51-P	Rat Uncoupling Protein 5 (UCP5) Control/blocking peptide
UCP52-A	Anti-Human Uncoupling Protein 5 (UCP5) IgG #2, aff pure
UCP52-P	Human Uncoupling Protein 5 (UCP5) Control/blocking peptide #2
UCP55-C	Recombinant purified Mouse Uncoupling Protein 5 (UCP5) (full length, His-tag) WB +ve control
UCP55-R-10	Recombinant purified Mouse Uncoupling Protein 5 (UCP5) (full length, His-tag)
UCP56-R-10	Recombinant purified Human Uncoupling Protein 5 (UCP5) (full length, His-tag)

UCP25-R-10-50 140121A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com