


# Antistatic Pour Boat, Small

Product ID W598

## Description

Polystyrene weigh boats designed to conveniently weigh liquids and statically-affected powders.

## Product Information

Autoclavability	Not autoclavable	
Capacity	20 mL	
Material	Polystyrene	
Product Number	W598	
Length	57 mm	
Width	44 mm	
Height	4 mm	

## Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	4823.20.9000

**PhytoTechnology Laboratories®**

14610 W 106th St. Lenexa, KS 66215

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com