

Product Specification Sheet

**Zinc Transporter 2 ((SLC39A2; Eti-1; 6A1)) Antibodies**

<b>Cat # ZIP21-P</b>	Mouse ZIP2 Control/Blocking Peptide	<b>SIZE:</b> 100 µg
<b>Cat # ZIP21-A</b>	Rabbit anti-mouse ZIP2 IgG (affinity pure)	<b>SIZE:</b> 100 µg

**Zinc** is an essential nutrient for all organisms because of the many important roles this metal plays. Movement of zinc into and out of cells and subcellular organelles is mediated by zinc transporter proteins. In many organisms, zinc uptake is mediated by members of the ZIP family of metal ion transporters. In mammals, the ZIP2, Zip2, Zip3, Zip4, Zip5, LIV-1 (Zip6), KE4 (Zip7), and BIGM103 (Zip8) proteins have been implicated in zinc uptake in a variety of cell and tissue types.

**ZIP2 ((SLC39A2; Eti-1; 6A1))** encodes a protein of 309 amino acids, includes an amino terminal signal peptide of 26 amino acids, and has the predicted topology of Zip proteins. The human or mouse *Zip2* overexpression in transfected K562 and human embryonic kidney (HEK293) cells, respectively, increased 65Zn uptake in these cells. *hZIP2*-induced Zn uptake was shown to be highly inhibited by Mn<sup>2+</sup> and Co<sup>2+</sup>, metals that did not substantially affect *hZIP1*-, *mZIP1*-, or *mZIP2*-induced Zn uptake. Zn uptake mediated by *hZIP2* appears to be stimulated by HCO<sub>3</sub><sup>-</sup>, non-ATP dependent, but sensitive to N-ethylmaleimide, a sulfhydryl-reactive agent.

**ZIP2** rat: 288aa; human: 309aa, mouse: 309aa – 33KDa; human Chromosome: 14q11.2. Expressed only in prostate and uterine epithelial cells.

**Source of Antigen, Antibodies**

<b>Antigen</b>	19- aa peptide of Mouse ZIP2 (Protein accession # <b>Q9NP94</b> ; ref. 1); designated as ZIP21-P control/blocking peptide conjugated to KLH; epitope location ~ N-terminus, Extracellular domain
<b>Antibody host/type</b>	Rabbit, Polyclonal IgG (Cat # ZIP21-A), purified over antigen-Agarose
<b>Secondary Ab</b>	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
<b>Negative Control Ab</b>	Non-immune rabbit IgG (Cat # 20009-1) to be used as –ve control for ELISA, WB, IHC etc.

**Form & Storage of Antibodies/Peptide Control**

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder** in PBS at 1mg/ml

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder** in PBS at 1 mg/ml.

**Storage**

**Short-term:** unopened, undiluted vials for less than a week at 4oC.

**Long-term:** at –20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at –20°C or below.

**Shipping:** 4°C for solutions and room temp for powder.

**Recommended Usage**

**Western Blotting:** 1-10 µg/ml; using affinity pure antibody (chemiluminescence technique).

**ELISA:** 1:100K; using 50-100 ng control peptide/well.

**Histochemistry & Immunofluorescence:** Not tested; we recommend the use of affinity purified antibody at 2-10 µg/ml.

**Specificity & Cross-reactivity**

Mouse ZIP21-P peptide sequence has no homology to mouse or rat ZIP2 protein. Antibody cross-reactivity in various species is not known. We recommend using antibody Cat # ZIP21-A against human ZIP2 protein. The control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity.

**General References:**

- 1) Yamaguchi S. (1995) Kokubyo Gakkai Zasshi 62:78-93

**List of related items, data sheets, and publications, using ADI antibodies is posted on the web site**

\*This product is for in vitro research use only.

**Related material available from ADI**

- Antibodies to human, mouse and rat ZIP1-7

ZIP21-A 71205A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)