



## Product Information Sheet

# B236 Buffer Solution, pH Reference Standard, pH 7.0 ± 0.02 at 25° C

### Properties

Form: Liquid  
Appearance: Clear, Yellow  
Application: pH Buffer Standard  
Solubility:  
Typical Working  
Concentration:  
Storage Temp: Room Temperature  
Storage Temp of  
Stock Solution:  
Other Notes:

### Application Notes

This is used for pH meter calibration

### References

Revised 3/2007

**PhytoTechnology Laboratories, LLC™**  
P.O. Box 12205 • Shawnee Mission, KS • 66282-2205™

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com