



Product Information Sheet

B237 **Buffer Solution, pH Reference** **Standard, pH 10.0 ± 0.02 at 25° C**

Properties

Form: Liquid
Appearance: Clear, Blue
Application: pH Buffer Standard
Solubility:
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes:

Application Notes

This is used for pH meter calibration

References

Revised 3/2007

PhytoTechnology Laboratories, LLC[™]
P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

India Contact:
Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com