



Product Information Sheet

M498 Malt Extract Agar

Properties

Form:	Powder
Appearance:	Off-White to Cream Powder
Application:	Bacterial Screening
Solubility:	Hot Water
Typical Working Concentration:	48.0 g/L
Storage Temp:	Room Temperature
Other Notes:	This product should be suspended in water and heated to boiling before dispensing and autoclaving.

Formula (mg/L)

Agar	15,000
Malt Extract	30,000
Peptone from Meat	3000

Application Notes

Typically used for the cultivation of yeasts and molds.

Revised Oct 2014

PhytoTechnology Laboratories, LLC

P.O. Box 12205; Shawnee Mission, KS 66282-2205

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com