




Product Information Sheet



S663
TRISODIUM CITRATE
DIHYDRATE

CAS: 6132-04-3
Formula: $C_6H_5Na_3O_7 \cdot 2H_2O$
Molecular Wt: 294.10

Properties

Form: Granular Powder
Appearance: White
Application: Generally used as a buffer
Solubility: Soluble in Water
Typical Working Concentration: N/A
Storage Temp: Room Temperature
Storage Temp of Stock Solution: Room Temperature

Application Notes

Commonly used for the preparation of molecular biology buffers.
Trisodium citrate dihydrate is a tribasic salt of citric acid.

PhytoTechnology Laboratories®
P.O. Box 12205; Shawnee Mission, KS 66282-2205

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com