

Product Information Sheet

A126 Antioxidant Mixture

Properties

Form:	Fine to Fluffy Powder
Appearance:	White
Application:	Antioxidant
Solubility:	Water
Typical Working Concentration:	0.25 g/L
Storage Temp:	Room Temperature
Storage Temp of Stock Solution:	2-6° C
Other Notes:	Plant Tissue Culture Tested

Formula (mg/L)

L-Ascorbic Acid	100
Citric Acid	150

Application Notes

This product typically is made up as a stock solution and an explant is dipped into it before transfer to fresh medium.

References

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

A126-Info

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400
Email - customerservice@lifetechindia.com | customerservice@atzlabs.com

Page 1 of 1