

Product Information Sheet

A143

L-Arginine, Free Base



Synonym: (S)-2-Amino-5-[(aminoiminomethyl)amino]pentanoic Acid

CAS: 74-79-3

Formula: C₆H₁₄N₄O₂

Molecular Wt: 174.20

Properties

Form: Powder

Appearance: White Crystalline or Granular

Application: Amino Acid

Solubility: Water

Typical Working

Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution:

Other Notes: Plant Tissue Culture Tested

May be Derived from Bacterial Fermentation of a Non-Animal Source
or Derived from an Avian Source

Application Notes

References

Merck 13, 785

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

A143-Info

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com

Page 1 of 1