

## Product Information Sheet

### B138 BM-1 Terrestrial Orchid Medium w/out Agar

#### Properties

Form:	Powder
Appearance:	White to Yellow Powder
Application:	Orchid Culture
Solubility:	Water
Typical Working Concentration:	21.22 g/L
Storage Temp:	2-6°
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insolubles may form.
Other Notes:	Same formulation as B141 without agar.

#### Formula

(mg/L)

Boric Acid	10
Cobalt Chloride•6H <sub>2</sub> O	0.025
Cupric Sulfate•5H <sub>2</sub> O	0.025
Na <sub>2</sub> EDTA•2H <sub>2</sub> O	37.25
Ferrous Sulfate•7H <sub>2</sub> O	27.85
Magnesium Sulfate, Anhydrous	100
Manganese Sulfate•H <sub>2</sub> O	25
Molybdc Acid (Sodium Salt) •2H <sub>2</sub> O	0.25
Potassium Phosphate, Monobasic	300
Zinc Sulfate•7H <sub>2</sub> O	10
D-Biotin	0.05

Casein, Enzymatic Hydrolysate	500
Folic Acid	0.5
L-Glutamine	100
Glycine (Free Base)	2
myo-Inositol	100
Nicotinic Acid (Free Acid)	5
Pyridoxine•HCl	0.5
Sucrose	20,000
Thiamine•HCl	0.5

#### Application Notes

Plant Tissue Culture Tested

Especially suited for terrestrial orchids, Paphiopedilum, and Phragmipedium.

Seed germination may be enhanced by the addition of 50 mL/L Coconut Water (Prod. No. C195).

Consult the Orchid Media Selection & Use Guide for additional information; this can be downloaded from the Technical Info page on the PhytoTechnology Labs web site.

Revised 12/2012

#### India Contact

### Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com