


Product Information Sheet

D217



6-(γ,γ -Dimethylallylamino)purine Solution (1 mg/mL)

Synonym: 2iP; N6-[2-Isopentyl]adenine
 CAS: 2365-40-4
 Formula: C₁₀H₁₃N₅
 Molecular Wt: 203.3

Properties

Form: Liquid
 Appearance: Colorless, Clear Liquid
 Application: Cytokinin
 Solubility: Miscible with Water
 Typical Working Concentration:
 Storage Temp: -20 to 0° C
 Storage Temp of Stock Solution:
 Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *PhytoTechnology Laboratories™* tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com