


Product Information Sheet

F323
Food Color, Red


Synonym: Erythrosine; FD&C Red #3

CAS: 16423-68-0

 Formula: C₂₀H₆I₄Na₂O₅

Molecular Wt: 879.86

Properties

Form: Powder

Appearance: Dark Red Powder

Application: Media Coloring

Solubility: Water

Typical Working Concentration: Recommended 10 milligrams to 250 milligrams per liter

Storage Temp: Room Temperature

Storage Temp of Stock Solution: Room Temperature

Other Notes: Plant Tissue Culture Tested

Application Notes

Once dissolved in water, this product produces a bright red color.

References

Merck 13, 3727

Revised 2/2016

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com