


Product Information Sheet

I409



Indole-3-Propionic Acid

Synonym: IPA
 CAS: 830-96-6
 Formula: C₁₁H₁₁NO₂
 Molecular Wt: 189.2

Properties

Form: Powder
 Appearance: Cream to Yellow Crystalline Powder
 Application: Auxin
 Solubility: KOH
 Typical Working Concentration:
 Storage Temp: 2 to 6° C
 Storage Temp of -20 to 0° C
 Stock Solution:
 Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *PhytoTechnology Laboratories™* tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com