

## Product Information Sheet

**M534**

### Malmgren Modified Terrestrial Orchid Medium Without Sucrose and Agar

#### Properties

Form:	Powder
Appearance:	Gray to Black
Application:	Orchid Culture
Solubility:	Water
Typical Working Concentration:	21.84 g/L
Storage Temp:	2 – 6°C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	Contains the macro- and micronutrients and organic constituents as described by Malmgren (1996). Does not contain Sucrose or Agar.

#### Formula

**(mg/L)**

Calcium Phosphate, Tribasic	75
Na <sub>2</sub> EDTA•2H <sub>2</sub> O	37.26
Ferrous Sulfate•7H <sub>2</sub> O	27.8
Magnesium Sulfate, Anhydrous	97.69
Manganese Sulfate•H <sub>2</sub> O	1.54
Potassium Phosphate, Monobasic	75
Activated Charcoal	1000
D-Biotin	0.05
Casein, Enzymatic Hydrolysate	400

Folic Acid	0.5
Glycine (Free Base)	2
myo-Inositol	100
Nicotinic Acid (Free Acid)	5
Pineapple Powder	20,000
Pyridoxine•HCl	5
Thiamine•HCl	10

#### Application Notes

Plant Tissue Culture Tested

Plant species: Terrestrial orchids including Paphiopedilum, Cypripedium, and Phragmipedium  
 This medium is a good seed sowing medium for many terrestrial orchid species. Originally developed for Cypripedium, this medium is also used for germination of Phragmipedium and Paphiopedilum. Germination may be enhanced by including 50 ml of Coconut Water (C195) per liter of medium.

#### References

Malmgren, S. 1996. In: C. Allen, Editor, North American Native Terrestrial Orchids. *Propagation and Production Conf. Proc.*, pp 63-71.

#### India Contact

### Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

M534-Info

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400  
 Email - customerservice@lifetechindia.com | customerservice@atzlabs.com

Page 1 of 1  
 Revised 5/2014