

Product Information Sheet

N611 Nutrient Broth

Properties

Form:	Powder
Appearance:	Cream to Tan
Application:	Microbiology
Solubility:	Water
Typical Working Concentration:	8.0 g/L
Storage Temp:	Room Temperature

Formula (mg/L)

Meat Extract	3,000
Peptone from Meat	5,000

Application Notes

This medium is used for the cultivation of a wide variety of microorganisms, most notably for nonfastidious organisms. The American Public Health Association (APHA) suggests nutrient broth as a standard culture medium for use in water testing (1917).

Nutrient broth is also a standard medium in many procedures for food testing as well as general microbiology applications. The meat extract and peptone in nutrient agar provide the nitrogen, carbon, vitamins and amino acids necessary for microorganism growth.

References

- American Public Health Association. 1917. Standard methods of water analysis, 3rd ed. American Public Health Association, Washington, D.C.
- Association of Official Analytical Chemists. 1995. Official methods of analysis of AOAC International, 16th ed. AOAC International, Arlington, VA.
- Bacteriological Analytical Manual. 1995. 8th ed. Association of Official Analytical Chemists. Gaithersburg, MD
- Eaton, A. D., L. S. Clesceri, and A. E. Greenberg (eds.). 1995. Standard methods for the examination of water and wastewater, 19th ed. American Public Health Association, Washington, D.C.
- Marshall, R. T. (ed.). 1993. Standard methods for the microbiological examination of dairy products, 16th ed. American Public Health Association, Washington, D.C.
- Vanderzant, C., and D. F. Splittstoesser (eds.). 1992. Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com