

Product Information Sheet

P6647**Phytotech Spathiphyllum
Medium****Properties**

Form:	Powder
Appearance:	White to Cream
Application:	Plant Tissue Culture
Solubility:	Partially Soluble in Water
Typical Working Concentration:	33.85 g/L
Storage Temp:	2 to 6 °C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	Proprietary formulation for the micropropagation of <i>Spathiphyllum</i> ; contains a gelling agent, carbohydrates, and plant growth regulators.

Application Notes

Plant Tissue Culture Tested

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com