


**Product Information Sheet****P725  
Polyethylene Glycol**

Synonym: PEG-1500; Carbowax

CAS: 25322-68-3

Formula:  $\text{H}(\text{C}_2\text{H}_4\text{O})_n\text{OH}$ 

Molecular Wt: 1300-1600

**Properties**

Form: Solid

Appearance: White Waxy Solid

Application: Molecular Biology

Solubility: Water

Typical Working

Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution:

Other Notes: Plant Tissue Culture Tested

**Application Notes****References**

Merck 13, 7651

Revised 3/2007

**India Contact****Life Technologies (India) Pvt Ltd.**

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com