

Product Information Sheet

KCl

P730

Potassium Chloride Solution 3.0 M

Synonym: Electrode Filling/ Soaking Solution
 CAS: 7447-40-7
 Formula: KCl
 Molecular Wt: 74.55

Properties

Form: Liquid
 Appearance: Colorless, Clear Solution
 Application: Electrode Filling/ Soaking
 Solubility: Miscible with Water
 Typical Working Concentration:
 Storage Temp: Room Temperature
 Storage Temp of Stock Solution: N/A
 Other Notes:

Formula (g/L)

Potassium Chloride	223.65
--------------------	--------

Application Notes

References

Merck 13, 7704

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com