


Product Information Sheet

P770
Pluronic® F-68



Synonym: Poloxamer 188

CAS: 9003-11-6

Formula: $(C_3H_6O \cdot C_2H_4O)_x$

Molecular Wt: 8400 Avg.

Properties

Form: Powder

Appearance: White Powder

Application: Surfactant

Solubility: Water

Typical Working

Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution:

Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck 13, 7645

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com