



PHYTO TECHNOLOGY LABORATORIES[®]

"Helping to Build a Better Tomorrow Through Plant Science"[™]

Glass Bead Sterilizer

ErgoSteri VT[™]



PhytoTechnology Laboratories is pleased to introduce their new ErgoSteri[™] glass bead sterilizer. This unique unit features a slanted bead well that is more ergonomically friendly than convention units that are vertical. Another important feature is a digital temperature controller that is adjustable between 110 and 350° C. Most units on the market today feature either factory pre-set temperatures that cannot be changed, or analog knobs that make it difficult to set an exact temperature. The only other digital model, the Steripot brand, has a very limited temp range (250 to 275° C).

The ErgoSteri[™] was designed with quality and accurate temperature maintenance in mind. The unit is manufactured and assembled in the USA, and the digital temperature controller allows for accurate temperature adjustments within a broad range. The construction of the unit allows for easy maintenance and care-free operation; however, *PhytoTech* offers a one-year warranty on parts & labor should a problem be encountered. And since all manufacturing and assembly are performed in the USA *PhytoTech* can offer a quick turnaround if any unit should require warranty work.

The ErgoSteri[™] heats fast, reading a stable 200° C in under 20 minutes. *PhytoTech's* unique design also prohibits glass beads from sliding into any other part of the unit other than the bead well making filling and emptying much easier.

Please visit our website for more information on the ErgoSteri[™] glass bead sterilizer.

Product Number	Product Description	Product Notes		Dimensions and Other Info		Pkg Size	Price (USD)*
S7510	ERGOSTERI VT [™] , 110 V Glass Bead Sterilizer	Temperature Read-Out: Digital Temperature Range: 110 to 350° C Manufactured in: USA Warranty: 1-year	110 to 350° C	Body Height:	20.1 cm (8.25 in)	1 ea	\$525.75
S7720	ERGOSTERI VT [™] , 220 V Glass Bead Sterilizer			Body Width:	14.5 cm (5.75 in)		
		Body Length (front to back):	23.5 cm (9.25 in)				
		Bead Well Diameter:	5 cm (2 in)				
		Bead Well Depth:	12.7 cm (5 in)				
S638	BEAD REPLACEMENTS Glass Bead Sterilizer Beads	Packaged specifically for the ErgoSteri bead well capacity.		Quantity:	450 g	1 ea	\$22.75

*prices are subject to change without notification.

Competitor	Model	Country of Origin	Temperature	Temperature Read-Out	Ergonomic Design
A	MODEL BS-1000	INDIA	PRE-SET	NONE	NO
B	STERI-350	SWITZERLAND	PRE-SET	ANALOG	NO
C	STERIPOT-S150D with Digital Controller	INDIA	ADJUSTABLE	DIGITAL	NO
D	LAB ASSOCIATES G3301	THE NETHERLANDS	ADJUSTABLE	NONE	NO
PhytoTech	ERGOSTERI VT[™]	USA	ADJUSTABLE	DIGITAL	YES