


Product Information Sheet**S807**
L-Serine

Synonym: (S)-2-Amino-3-hydroxypropanoic Acid

CAS: 56-45-1

Formula: C₃H₇O₃N

Molecular Wt: 105.09

Properties

Form: Powder

Appearance: White Crystalline

Application: Amino Acid

Solubility: Water

Typical Working

Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution:

Other Notes: Plant Tissue Culture Tested

Derived from Bacterial Fermentation of a Non-animal Source

Application Notes**References**

Merck 13, 8534

Revised 3/2007

India Contact**Life Technologies (India) Pvt Ltd.**

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com