

**Product Information Sheet****T7900**  
**Tryptic Soy Broth I**

Synonym: TSB

**Properties**

Form: Powder  
Appearance: Cream to Tan Powder  
Application: Bacterial Screening, Microbiology & Phytopathology  
Solubility: Soluble in Water  
Typical Working Concentration: 25.00 g/L  
Storage Temp: Room Temperature  
Storage Temp of Stock Solution: Preparation of concentrated solutions is not recommended as insoluble precipitates may form.

**Formula****(g/L)**

Formula	(g/L)
Tryptone	15
Peptone from Soy	3
Sodium Chloride	5

**Application Notes**

This is a general purpose medium for the cultivation of fastidious microorganisms.

Please Note: It is the sole responsibility of the purchaser to determine the appropriateness of this product.

**References**

1. Difco Manual, 1984. Pg. 1027-1028.

**India Contact****Life Technologies (India) Pvt Ltd,**

T7900-Info

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)  
Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400  
Email - customerservice@lifetechindia.com | customerservice@atzlabs.com

Page 1 of 1  
Oct 2014