

Product Information Sheet

V891

Vacin & Went Modified Orchid Medium

Properties

Form:	Powder
Appearance:	White to Off-White powder
Application:	Orchid Culture
Solubility:	Water
Typical Working Concentration:	21.67 g/L
Storage Temp:	2 – 6° C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	Contains the macro- and micronutrients as described by Vacin and Went (1949) modified with an equivalent iron molar concentration of ferrous sulfate in place of ferric tartrate. Contains Sucrose

Formula

Ammonium Sulfate	500
Calcium Phosphate, Tribasic	200
Na ₂ EDTA•2H ₂ O	37.26
Ferrous Sulfate•7H ₂ O	27.8
Magnesium Sulfate, Anhydrous	122.1
Manganese Sulfate•H ₂ O	5.6
Potassium Nitrate	525
Potassium Phosphate, Monobasic	250
Sucrose	20,000
Thiamine•HCl	0.4

Application Notes

Plant Tissue Culture Tested

References

Vacin, EF and EW Went. 1949. Bot. Gaz. 110: 605-613.

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com