


Product Information Sheet

Z875



Zeatin Riboside Solution (1 mg/mL), *trans* Isomer

Synonym: N/A
CAS: 6025-53-2
Formula: C₁₅H₂₁N₅O₅
Molecular Wt: 351.4

Properties

Form: Liquid
Appearance: Colorless, Clear Liquid
Application: Cytokinin
Solubility: Miscible with Water
Typical Working Concentration:
Storage Temp: -20 to 0° C
Storage Temp of Stock Solution:
Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *PhytoTechnology Laboratories™* tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References

India Contact

Life Technologies (India) Pvt Ltd,

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com